

SiLA Rapid Integration[®]

Introducing SiLA 2



Who is he?



Oskari Vinko



Who is SiLA?

- Non-profit organisation founded in 2008
- Global footprint with 35 member companies
- 1 600 + individual members



SiLA Mission

- Create international, open connectivity standard in lab automation
- Improve data integrity in labs
- Accessible, straightforward, stable, cost effective, usable and extensible

You spoke

- We restricted you to predefined device types
- Connecting a new device wasn't easy
- SiLA 1 was based on XML/SOAP
- Standardisation process was slow
- It was hard to get started!

We listened

SiLA 2

SiLA 2 Device

Provides features



Plate reading
Plate shaking
Temperature control

SiLA 2 Device

Wireless or wired



SiLA 2 Device

Discover automagically!



SiLA 2 Device


Discover automagically!



Plate reading
Plate shaking
Temperature control

SiLA 2 Device

Peer to peer



SiLA 2 Device

Subscribe to properties




37°C




SiLA 2 Device

Subscribe to properties




SiLA 2 Device

Subscribe to properties



SiLA 2 Device

Subscribe to properties



SiLA 2

Out late 2018

Release candidate accessible to members in spring 2018

Vendor is the designer

- Use approved features or design your own
- Connect other devices for synergistic features
 - Barcode reader + LIMS + Liquid handler
- Benefit from cross-platform and cross-language libraries to make driver development a breeze!
 - The libraries will be open source
 - Come to our booth, we'll show you personally how!
- Increase customer value via better usability and compatibility

User has the power

- *Ad hoc* automation: in the near future, anyone can configure simple automation themselves
- Break free from closed vendor ecosystems
 - Re-use the SiLA 2 devices in new automation projects
- Quicker and cheaper integration
 - Get your lab connected to your IT systems
- Benefit from open source SiLA 2 driver repository
 - Share the drivers you paid for
 - Download SiLA 2 drivers for free
 - Buy converters from member companies
- Try it out at our booth!


New integration opportunities

- Use a standardized device network designed **with/having?** easy integration in mind
 - Save effort on device drivers, feature design, device interaction, real-time communication, discovery etc.
- Design device agnostic software
 - User empowering solutions that let them decide
 - Test the SiLA Browser in action and get inspired!
- Bring convenience to the user

Ecosystem for everyone!

- Several applications can use the SiLA peer to peer network simultaneously
- Exciting opportunity for startups:
 - Access ready SiLA ecosystem and infrastructure
- No license fees, no royalties, no membership
 - It's simply free!
- SiLA certification is optional, proof of quality

Ecosystem for everyone!



SiLA 2 Browser demo



- Here we will add:
 - A demo video or screenshots of the open source SiLA Browser (developed by UniteLabs)
 - Same demo is available at SiLA booth
 - See SiLA 2 capabilities in action (discovery, device agnostic commanding etc.)

How to get started?

- Join SiLA now, or wait until end of 2018
 - XX % discount during SLAS
- Participate in training programs provided by SiLA
- Wait for our high quality web resources
 - SiLA 2 specification sheet
 - Getting started guide
 - SiLA browser and other software components
 - Open source libraries with crystal clear tutorials

How to join?

- Individuals: register free on our website
- Organisations: fill out the membership form at
 - <http://www.sila-standard.com/join-us/>
 - Membership fee starts from 3500 EUR / year
- Benefits:
 - Network with end users, vendors, integrators and consultants
 - Get support and access expertise
 - Visibility and collaboration opportunities
 - Get an early access to SiLA 2!

How to contribute?

- For SiLA members:
 - Participate to the SiLA 2 work group
 - daniel.juchli@silastandard.org
 - Specification, documentation, reference implementation
 - Team up with other members to upgrade your lab/product to SiLA 2
- For non-members:
 - Wait for public release
 - Share your SiLA 2 drivers on the open driver repository
 - Contribute to the SiLA 2 libraries on GitHub

SiLA Rapid Integration®

Thank you!

SiLA 2 will be out late 2018,
release candidate in spring 2018

<http://www.sila-standard.org/>